


Potato Virus X ELISA Kit Cat# POTX-096

The Potato Virus Y ELISA Kit is a highly sensitive sandwich ELISA for the measurement of Potato Virus Y in plant extracts



**Potato Virus X ELISA Kit features**

- Anti-PVX antibody pre-coated, ready-to-use 96-well breakable strip plate, suitable for multiple runs over 6-12 months
- Recombinant PVX Coat protein standard
- **Assay length:** 3 hours & 45 minutes. 4 incubation steps at room temperature
- **Dynamic range:** 156-10,000 pg/ml
- **Sensitivity** ~40 pg/ml
- **Storage:** 2-8°C (whole kit)
- **Shelf life:** 6-12 months

Contains all necessary reagents *For in vitro research use only*

**Assay Procedure:** Allow all reagents to reach room temperature. Arrange and label required number of strips.

- Step 1.** Pipette 100 ul of appropriately diluted samples and calibrators into wells and incubate for 2 hours at room temperature.
- Step 2.** Wash the wells 3X with 300 ul of wash buffer for each well
- Step 3.** Add 100 ul of Biotin conjugated detection antibody to each well and incubate for 1 hour at room temperature
- Step 4.** Wash the wells 3X with 300 ul of wash buffer for each well
- Step 5.** Add 100 ul of Streptavidin HRP conjugate to each well and incubate for 30 minutes at room temperature
- Step 6.** Wash the wells 3X with 300 ul of wash buffer for each well
- Step 7.** Add 100 ul of TMB Substrate solution to all wells, mix gently, and incubate at room temperature for 15 minutes.
- Step 8.** Pipette 100 ul of stop solution into each well and mix gently. Measure at 450 nm w/ 630 nm as a reference filter if available.

**Performance Characteristics**

**Sensitivity:** ~40 pg/ml  
**Average recovery:** 100 ±15%  
**Average linearity:** 100 ±15%  
**Precision:** Intra-assay: <10%      Inter-assay: <10%

**Minimum recommended dilution**  
**Plant extracts: 20-fold**

**Note:** Minimum recommended dilution represents the dilution which is needed to eliminate matrix interference effects and obtain optimal recovery. All samples must be diluted to at least the minimum recommended ratio. Samples may be further diluted if the sample values fall within the standard curve, if sample volume is to be preserved, or if the sample value is above the highest OD on the standard curve

**General Information**

Potato virus X (PVX) is a pathogenic plant virus belonging to the family Alphaflexiviridae. PVX is primarily found in potatoes and is only transmitted mechanically. The virus causes mild to no symptoms in most potatoes. When the potato is infected with both Potato Virus Y and Potato Virus X, severe symptoms can be caused in potatoes. The virus is a single-stranded positive-sense RNA virus that is wrapped in approximately 1300 units of a single coat protein